



Plaster for bandages on the basis of synthetic alpha calcium sulphate hemihydrate

Charakteristics

Sieve residue on 500 µm	max. 0,1 %
Crystal water	max. 6,2 %
Degree of whiteness (powder)	min. 94 %
Final setting time (at mixing ratio powder:water 100:50 pbw.)	3,5 - 5 min
Tensile strength (moist) (at mixing ratio powder:water 100:67 pbw.)	min. 4 N/mm ²

Packing

Paperbags with foil liner 50 kgs net

Storage

At least 12 month in closed, vapour-proof packages.

The recommendations are given to the best of our knowledge after careful control. We grant the quality of our products according to our specification. Any further liability cannot be accepted since the proper application of our products is outside of our control.